

# LSM14A Antibody



PACO10300

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

LSM14A

**Uniprot**

Q8ND56

**Synonyms:**

LSM14A, SCD6 homolog A (S. cerevisiae); LSM14A; C19orf13; DKFZp434D1335; FAM61A; RAP55 ;

**Immunogen:**

Human LSM14A.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A