

LSM5 Antibody



PACO10303

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

LSM5

Uniprot

Q9Y4Y9

Synonyms:

LSM5 homolog, U6 small nuclear RNA associated (S. cerevisiae); LSM5; FLJ12710; YER146W ;

Immunogen:

Human LSM5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A