

LY96 Antibody



PACO10314

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

LY96

Uniprot

Q9Y6Y9

Synonyms:

lymphocyte antigen 96; LY96; MD-2; MD2 ;

Immunogen:

Human LY96.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A