

LYZL6 Antibody



PACO10321

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

LYZL6

Uniprot

O75951

Synonyms:

lysozyme-like 6; LYZL6;1700023H08Rik; LYC1; PRO1485; TKAL754 ;

Immunogen:

Human LYZL6.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A