

LZTS1 Antibody



PACO10324

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:**Protein Background:****Gene ID:**

LZTS1

Uniprot

Q9Y250

Synonyms:

leucine zipper, putative tumor suppressor 1; LZTS1; F37; FEZ1 ;

Immunogen:

Human LZTS1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A