

MAD2L1BP Antibody



PACO10328

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

MAD2L1BP

Uniprot

Q15013

Synonyms:

MAD2L1 binding protein; MAD2L1BP; CMT2; KIAA0110; MGC11282; RP1-261G23.6 ;

Immunogen:

Human MAD2L1BP.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A