

MAEA Antibody



PACO10330

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

MAEA

Uniprot

Q7L5Y9

Synonyms:

macrophage erythroblast attacher; MAEA; EMLP; EMP; HLC-10; PIG5 ;

Immunogen:

Human MAEA.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A