

MAGEA8 Antibody



PACO10338

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MAGEA8

Uniprot

P43361

Synonyms:

melanoma antigen family A, 8; MAGEA8; CT1.8; MAGE8; MGC2182 ;

Immunogen:

Human MAGEA8.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A