

MAGEC2 Antibody



PACO10342

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

MAGEC2

Uniprot

Q9UBF1

Synonyms:

melanoma antigen family C, 2; MAGEC2; CT10; HCA587; MAGEE1; MGC13377 ;

Immunogen:

Human MAGEC2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A