

MAGI1 Antibody



PACO10347

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MAGI1

Uniprot

Q96QZ7

Synonyms:

membrane associated guanylate kinase, WW and PDZ domain containing 1; MAGI1; AIP3; BAIAP1; BAP1; MAGI-1; TNRC19; WWP3 ;

Immunogen:

Human MAGI1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A