

MAK16 Antibody



PACO10351

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

MAK16

Uniprot

Q9BXY0

Synonyms:

MAK16 homolog (S. cerevisiae); MAK16; MAK16L; RBM13 ;

Immunogen:

Human MAK16.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A