

MAN1A2 Antibody



PACO10353

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

MAN1A2

Uniprot

O60476

Synonyms:

mannosidase, alpha, class 1A, member 2; MAN1A2; MAN1B ;

Immunogen:

Human MAN1A2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A