

MAP1LC3B Antibody



PACO10359

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

MAP1LC3B

Uniprot

Q9GZQ8

Synonyms:

microtubule-associated protein 1 light chain 3 beta; MAP1LC3B; LC3B; MAP1A/1BLC3 ;

Immunogen:

Human MAP1LC3B.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A