

MAP4 Antibody



PACO10371

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

MAP4

Uniprot

P27816

Synonyms:

microtubule-associated protein 4; MAP4; DKFZp779A1753; MGC8617 ;

Immunogen:

Human MAP4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A