

MAPK11 Antibody



PACO10378

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

MAPK11

Uniprot

Q15759

Synonyms:

mitogen-activated protein kinase 11; MAPK11; P38B; P38BETA2; PRKM11; SAPK2; SAPK2B; p38-2; p38beta ;

Immunogen:

Human MAPK11.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A