

MAPK6 Antibody



PACO10380

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

MAPK6

Uniprot

Q16659

Synonyms:

mitogen-activated protein kinase 6; MAPK6; DKFZp686F03189; ERK3; HsT17250; PRKM6; p97MAPK ;

Immunogen:

Human MAPK6.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A