

LAMTOR3 Antibody



PACO10387

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

LAMTOR3

Uniprot

Q9UHA4

Synonyms:

MAPK scaffold protein 1; MAPKSP1; MAP2K1IP1; MAPBP; MP1 ;

Immunogen:

Human MAPKSP1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A