

MAPRE2 Antibody



PACO10389

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

MAPRE2

Uniprot

Q15555

Synonyms:

microtubule-associated protein, RP/EB family, member 2; MAPRE2; EB1; EB2; RP1 ;

Immunogen:

Human MAPRE2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A