

MARCKS Antibody



PACO10391

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

MARCKS

Uniprot

P29966

Synonyms:

myristoylated alanine rich protein kinase C substrate; Marcks; Macs; PKCSL ;

Immunogen:

Human Marcks.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A