

MAST1 Antibody



PACO10397

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

MAST1

Uniprot

Q9Y2H9

Synonyms:

microtubule associated serine/threonine kinase 1; MAST1; KIAA0973; SAST; SAST170 ;

Immunogen:

Human MAST1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A