

MBD2 Antibody



PACO10405

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MBD2

Uniprot

Q9UBB5

Synonyms:

methyl-CpG binding domain protein 2; MBD2; DKFZp586O0821; DMTase; NY-CO-41 ;

Immunogen:

Human MBD2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A