

MBD5 Antibody



PACO10408

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

MBD5

Uniprot

Q9P267

Synonyms:

methyl-CpG binding domain protein 5; MBD5; FLJ11113; FLJ30517; KIAA1461; MRD1 ;

Immunogen:

Human MBD5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A