

MCM2 Antibody



PACO10419

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

MCM2

Uniprot

P49736

Synonyms:

minichromosome maintenance complex component 2; MCM2; BM28; CCNL1; CDCL1; D3S3194; KIAA0030; MGC10606; MITOTIN; cdc19 ;

Immunogen:

Human MCM2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A