

# MCM3 Antibody



PACO10420

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

MCM3

**Uniprot**

P25205

**Synonyms:**

minichromosome maintenance complex component 3; MCM3; HCC5; MGC1157; P1-MCM3; P1. h; RLFB ;

**Immunogen:**

Human MCM3.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A