

# MCM4 Antibody



PACO10421

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

MCM4

**Uniprot**

P33991

**Synonyms:**

minichromosome maintenance complex component 4; MCM4; CDC21; CDC54; MGC33310; P1-CDC21; hCdc21 ;

**Immunogen:**

Human MCM4.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A