

MCM8 Antibody



PACO10425

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MCM8

Uniprot

Q9UJA3

Synonyms:

minichromosome maintenance complex component 8; MCM8; C20orf154; MGC119522; MGC119523; MGC12866; MGC4816; REC; dJ967N21.5 ;

Immunogen:

Human MCM8.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A