

ME1 Antibody



PACO10433

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

ME1

Uniprot

P48163

Synonyms:

malic enzyme 1, NADP(+)-dependent, cytosolic; ME1; HUMNDME; MES ;

Immunogen:

Human ME1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A