

MED12 Antibody



PACO10437

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MED12

Uniprot

Q93074

Synonyms:

mediator complex subunit 12; MED12; CAGH45; FGS1; HOPA; KIAA0192; OKS; OPA1; TNRC11; TRAP230 ;

Immunogen:

Human MED12.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A