

MED12L Antibody



PACO10438

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MED12L

Uniprot

Q86YW9

Synonyms:

mediator complex subunit 12-like; MED12L; FLJ46288; KIAA1635; NOPAR; TNRC11L; TRALP; TRALPUSH ;

Immunogen:

Human MED12L

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A