

MED17 Antibody



PACO10440

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MED17

Uniprot

Q9NVC6

Synonyms:

mediator complex subunit 17; MED17; CRSP6; CRSP77; DRIP80; FLJ10812; TRAP80 ;

Immunogen:

Human MED17.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A