

# MED4 Antibody



PACO10447

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

MED4

**Uniprot**

Q9NPJ6

**Synonyms:**

mediator complex subunit 4; MED4; DRIP36; FLJ10956; HSPC126; RP11-90M2.2; TRAP36; VDRIP ;

**Immunogen:**

Human MED4.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A