

METAP2 Antibody



PACO10460

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

METAP2

Uniprot

P50579

Synonyms:

methionyl aminopeptidase 2; METAP2; MAP2; MNPEP; p67; p67eIF2 ;

Immunogen:

Human METAP2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A