

MFSD7 Antibody



PACO10475

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

SLC49A3

Uniprot

Q6UXD7

Synonyms:

major facilitator superfamily domain containing 7; MFSD7; FLJ22269; LP2561 ;

Immunogen:

Human MFSD7.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A