

# MGEA5 Antibody



PACO10481

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:****Protein Background:****Gene ID:**

MGEA5

**Uniprot**

O60502

**Synonyms:**

meningioma expressed antigen 5 (hyaluronidase); MGEA5; FLJ11229; FLJ23355; KIAA0679; MEA5; NCOAT; OGA ; O-GlcNAcase

**Immunogen:**

Human MGEA5.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A