

MGMT Antibody



PACO10483

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

MGMT

Uniprot

P16455

Synonyms:

O-6-methylguanine-DNA methyltransferase;

Immunogen:

Human MGMT.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A