

MICAL2 Antibody



PACO10489

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

MICAL2

Uniprot

O94851

Synonyms:

microtubule associated monooxygenase, calponin and LIM domain containing 2; MICAL2; DKFZp686H03148; DKFZp686H2469; KIAA0750; MICAL2PV1; MICAL2PV2 ;

Immunogen:

Human MICAL2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A