

MIER2 Antibody



PACO10492

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MIER2

Uniprot

Q8N344

Synonyms:

mesoderm induction early response 1, family member 2; MIER2; KIAA1193 ;

Immunogen:

Human MIER2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A