

MIER3 Antibody



PACO10493

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MIER3

Uniprot

Q7Z3K6

Synonyms:

mesoderm induction early response 1, family member 3; MIER3; DKFZp686L09111; DKFZp781G11119; DKFZp781I11119; FLJ35954 ;

Immunogen:

Human MIER3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A