

MINK1 Antibody



PACO10496

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MINK1

Uniprot

Q8N4C8

Synonyms:

misshapen-like kinase 1 (zebrafish); MINK1; B55; MAP4K6; MGC21111; MINK; YSK2; ZC3; hMINK; hMINKbeta ;

Immunogen:

Human MINK1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A