

MKLN1 Antibody



PACO10503

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MKLN1

Uniprot

Q9UL63

Synonyms:

muskelin 1, intracellular mediator containing kelch motifs; MKLN1; FLJ11162; TWA2 ;

Immunogen:

Human MKLN1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A