

# MLLT4 Antibody



PACO10512

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

AFDN

**Uniprot**

P55196

**Synonyms:**

myeloid/lymphoid or mixed-lineage leukemia (trithorax homolog, Drosophila);  
translocated to, 4; MLLT4; AF-6; AF6; AFADIN; FLJ34371; RP3-431P23.3 ;

**Immunogen:**

Human MLLT4.

**Storage:**

0.1M NaHCO<sub>3</sub>, 0.1M Glycine, 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C; C,  
Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A