

MMP19 Antibody



PACO10524

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MMP19

Uniprot

Q99542

Synonyms:

matrix metalloproteinase 19; MMP19; MMP18; RASI-1 ;

Immunogen:

Human MMP19.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C; Avoid freeze / thaw cycles.

Product Images

N/A

N/A