

# MMP28 Antibody



PACO10526

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

MMP28

**Uniprot**

Q9H239

**Synonyms:**

matrix metalloproteinase 28; MMP28; EPILYSIN; MM28; MMP-28; MMP25 ;

**Immunogen:**

Human MMP28.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A