

# MMP7 Antibody



PACO10528

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

MMP7

**Uniprot**

P09237

**Synonyms:**

matrix metalloproteinase 7 (matrilysin, uterine); MMP7; MMP-7; MPSL1; PUMP-1 ;

**Immunogen:**

Human MMP7.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A