

MNDA Antibody



PACO10534

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MNDA

Uniprot

P41218

Synonyms:

myeloid cell nuclear differentiation antigen; MNDA; PYHIN3 ;

Immunogen:

Human MNDA.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A