

MOB1B Antibody



PACO10535

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

MOB1B

Uniprot

Q7L9L4

Synonyms:

MOB1, Mps One Binder kinase activator-like 1A (yeast); MOBKL1A; MATS2; MGC33910; MOB4A; Mob1B ;

Immunogen:

Human MOBKL1A.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A