

MOB1A Antibody



PACO10536

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MOB1A

Uniprot

Q9H8S9

Synonyms:

MOB1, Mps One Binder kinase activator-like 1B (yeast); MOBKL1B; C2orf6; FLJ10788; FLJ11595; MATS1; MOB1; MOBK1B; Mob4B ;

Immunogen:

Human MOBKL1B.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A