

MOB3A Antibody



PACO10537

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

MOB3A

Uniprot

Q96BX8

Synonyms:

MOB1, Mps One Binder kinase activator-like 2A (yeast); MOBKL2A; moblak ;

Immunogen:

Human MOBKL2A.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A