

MOCS2 Antibody



PACO10541

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MOCS2

Uniprot

O96033

Synonyms:

molybdenum cofactor synthesis 2; MOCS2; MCBPE; MOCO1; MOCS2A; MOCS2B; MPTS
;

Immunogen:

Human MOCS2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A