

MOG Antibody



PACO10542

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MOG

Uniprot

Q16653

Synonyms:

myelin oligodendrocyte glycoprotein; MOG; MGC26137; MOGIG-2 ;

Immunogen:

Human MOG.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A